



## GivMohr Sling Sizing Guide

Size	Height Range	Weight Range
Pediatric		<b>45-70 lb</b> (20.4-31.8 kg)
Extra-Small	<b>less than 5'0"</b> (152 cm)	<b>80-100 lb</b> (36.3-45.4 kg)
Small	<b>5'0"-5'4"</b> (152-163 cm)	<b>100-140 lb</b> (45.4-63.5 kg)
Medium	<b>5'3"-5'10"</b> (160-178 cm)	<b>135-185 lb</b> (61.2-83.9 kg)
Large	<b>5'10"-6'4"</b> (178-193 cm)	<b>180-235 lb</b> (81.6-106.6 kg)
Extra-Large	<b>over 6'4"</b> (193 cm)	<b>Over 235 lb</b> (106.6 kg)

### NOTES:

**GivMohr Sling** (1 arm) can be worn on either right or left arm.

**GivMohr Bilateral Sling** (2 arm) is worn on both arms simultaneously.

**Sizing Tip 1:** If "Height Range" and "Weight Range" yield different results, look at where the user carries their weight to assist in determining the proper size sling to use. If the weight is carried higher on the body (chest/shoulders), give priority to the "Weight Range" result. If the weight is carried lower on the body (stomach/hips), give priority to the "Height Range" result.

**Sizing Tip 2:** Custom sizing is available by mixing different sizes of the shoulder harness and arm strap. If you wish to explore this option, please contact us by email **BEFORE ORDERING** and provide the height and weight of the user.

