



GivMohr™ SLING

Flaccid Upper Extremity Positioning Device

Common Modifications to the GivMohr™ Sling

- [Securing the 4 Tab Buckle for optimal fit](#)
- [Adding wrist splint for proper wrist position](#)
- [Adjusting handle for arthritic thumbs](#)
- [Decreasing handle volume to accommodate distal motor return](#)
- [Mix sling sizes to accommodate different body types and arm lengths \(Bilateral\)](#)

Securing the 4 Tab Buckle for optimal fit

Step #1

Fit sling to user

Step #2

Remove sling carefully so 4 Tab Buckle does not slide

Step #3

Sew webbing straps together where they cross

[back to top](#)



Adding wrist splint for proper wrist position



[back to top](#)

Adjusting handle for arthritic thumbs

Step #1

Place sling in warm oven (lowest setting) until handle is soft (5-10 minutes)

Step #2

Pinch handle and cool in COLD water

[back to top](#)



Decreasing handle volume to accommodate distal motor return

Step #1

Cut plastic handle along webbing edge

Step #2

Wrap handle with mole skin to reduce handle size and pad rough edges.



[back to top](#)

Mix sling sizes to accommodate different body types

"Mix and Match" sizes of shoulder harness and elastic arm piece to accommodate unique sizes and/or different arm lengths (bilateral) of individuals

[back to top](#)



GivMohr Corporation
10000 Rio Grande BLVD N.W.
Albuquerque, N.M. 87114
Phone (505) 292-1144
Fax (505) 323-9526

© 2005-2009 GivMohr Corporation
GivMohr is a trademark of GivMohr Corporation